


PRODUCE PEDIGREES ON DEMAND

Progeny Anywhere is the pedigree drawing technology used in our software. We simply removed the control and license it directly to organizations or software vendors with programming resources. The control comes in Java to meet the needs of your current computing environment. Feed your current data into the Progeny Anywhere control via standard XML format. The control formats the data and generates a pedigree, providing the same features and functionality available within our software.

Displaying pedigrees in an easily readable and presentable diagram is a complex programming feat. Licensing our technology can give you the edge you need when it comes to displaying pedigrees.



YOUR DATA YOUR WAY

The data behind the pedigree can come from anywhere. The resulting pedigree can be customized the way you want. Utilize pedigree icon symbols and icon subtext to visualize your data.

Features	Benefits
Embedded Pedigree Technology	Add pedigree technology to your existing applications without having to write any code.
Pedigrees on Demand	Display pedigrees directly from your database with the click of a button.
Edit the Pedigree	Make changes to the pedigree and save those changes back to your database.
XML Interface	Exchange data via standard XML format; documentation explains exactly how to setup your programming calls.
Cost Effective	Save programming time by implementing a proven pedigree technology that has been in use commercially over the last 14 years.

PRICING

The Progeny Anywhere development license is offered on an academic or commercial basis. There is a cost for the development license and a per seat deployment charge. Please contact us for a price quote to meet your specific environment.

SYSTEM REQUIREMENTS

- Microsoft® Windows® 7, Vista, XP, 2000 (SP4) or Apple Mac OS 10.3 or later
- Intel® Pentium® III 500 MHz or compatible or Apple eMac 1.0GHz Motorola PowerPC G4 processor
- 64 MB RAM (recommended minimum); 256 MB RAM (recommended)
- 10 - 12 MB free hard disk space minimum

For more information, call:
800-PROGENY Toll Free
574-968-0822 Phone
574-968-0818 Fax

©Progeny Software, LLC. All Rights Reserved.
130 S. Main Street, Suite 420, South Bend, IN 46601

info@progenygenetics.com
www.progenygenetics.com