


CLINICAL WEB



CLINICAL WEB CONFIGURATION OPTION

Progeny Clinical Web offers you the flexibility to expand your system - connect from anywhere using any web browser on Windows or Mac. Clinical Web is an ideal cost-effective option for collaborative users, multi-platform users or users who are not tied to a single machine. Reach more users, reduce costs and simplify your updates - contact a sales rep to discuss how you can get started with Clinical Web today.



DATA SECURITY AND CONFIGURATION

Two different servers house the web server and database server separately. Web browsers will connect through an outside firewall to the Progeny Web Server which will securely direct traffic in and out of the database. The Progeny database is located behind an additional firewall, providing another security layer to prevent attacks on the web server from directly compromising the database.

CONNECT ANYTIME FROM ANY BROWSER

Run Progeny from any browser (IE, Firefox, Safari, Chrome) without compromising security and maintaining your HIPAA compliancy. Perform all the major functions as in a client installation (except for some administrative tasks).

SYSTEM REQUIREMENTS

The Progeny Web Server has the following requirements:

- Windows Server 2003/2008 Operating System
- 1GB RAM (2GB recommended)
- 100 MB Hard Drive Space
- Tomcat 6 (Servlet container/Application Server)
- Java 6

At least one Client installation is needed to perform the following administrative tasks:

- Field Creation
- Security Configuration
- Audit Setup
- Form Design
- Data Imports

For more information, call:
800-PROGENY Toll Free
574-968-0822 Phone
574-968-0818 Fax

©Progeny Software, LLC. All Rights Reserved.
130 S. Main Street, Suite 420, South Bend, IN 46601

info@progenygenetics.com
www.progenygenetics.com

