





## SECURITY AND COMPLIANCY

Features	Benefits
User and Class Based Security	User and class based security applies restrictions on what each user can perform or view in the sample lifecycle, and can also be set at the field, folder and database levels.
Set Restrictions on Sample Use	Define custom restrictions on how each sample is used for compliancy.
Complete Chain of Custody Control	Control when chain of custody is tracked for each custom-defined action performed on a sample throughout that sample's life cycle in the lab. Tracks what, who, when and what order each action was performed on a sample.
Custom Chain of Custody Reporting	Run reports on your chain of custody data for any time period. Reports can be sorted and filtered based on user, action, workflow and/or sample.
Compliance	Meets compliancy regulations, including GLP, 21 CFR Part 11, and HIPAA.
HIPAA Compliant Security Access Levels	Allow or prevent access to folders, fields, and samples. Auditing options allow for tracking passwords, logins, previous entry, new entry and user info with a date/time stamp.

## CORE DATABASE FUNCTIONALITY

Features	Benefits
Integration with Clinical and Lab	Progeny LIMS can be integrated seamlessly with Progeny Clinical (pedigree and clinical data management) and/or Progeny Lab (genotype data management) within a single application and using the same familiar interface.
Client and/or Browser Based Installation	Robust Windows based multi-user client for centralized database access. Platform independent browser option allows access from any location.
Create Unlimited Fields	Create any field you desire to track your data. Field types include text, date, number, yes/no, drop down, image, tables, computed, summary and hyperlink.
Create and Design Custom Data Entry Forms	Drag and drop to design custom data entry forms using your created fields. Setup multiple tabs, headers, questions, fields, hyperlinks, and more.
Enter, Search and Query Data	Search, query, view, sort, modify or enter data quickly in the Spreadsheet or through data entry forms.
Create Custom Reports	Custom reports can be generated and saved in just a few mouse clicks. Easily summarize your sample data based on any criteria or query.
Import and Export Data Easily	Import data from standard text files with user-defined formats, and export data to any text format. Customize and save export formats.
Data Validation	Set rules for data entry on a per-field basis to minimize errors and to ensure data integrity.

## SYSTEM REQUIREMENTS

**Client PC:** Windows 2000/XP/Vista; 1GB RAM recommended; 150MB hard drive space

**Server:** Windows 2000/XP/2003/2008, Linux, or Novell NetWare; 2GB RAM recommended, 200MB hard drive space required

**Web Browser:** Internet Explorer, Firefox, Safari, etc.

**Server:** Windows 2000/XP/2003/2008; 2GB RAM recommended, 200MB hard drive space required

**For more information, call:**  
800-PROGENY Toll Free  
574-968-0822 Phone  
+44 1379 854250 Europe  
574-968-0818 Fax

info@progenygenetics.com  
www.progenygenetics.com

©Progeny Software, LLC. All Rights Reserved.  
130 S. Main Street, Suite 420, South Bend, IN 46601 USA

